English	Dental Implant Surgical Instruments	2-3
Español	Instrumentos quirúrgicos para implantes dentales	4-5
Italiano	Strumenti chirurgici per impianti dentali	6-7
Português	Instrumentos cirúrgicos para implantes dentários	8-9
Deutsch	Chirurgische Instrumente für Zahnimplantate	10-11
Français	Instruments chirurgicaux pour implants dentaires	12-13
Nederlands	Chirurgische instrumenten voor dentale implantaten	14-15
Polski	Instrumenty chirurgiczne do implantów dentystycznych	16-17

Bicon, Llc
501 Arborway • Boston, MA 02130 USA

Bicon Europe, Ltd.
Unit 4 Ballycummin Village • Ballycummin • Limerick, Ireland

Description

Bicon Dental Implant Surgical Instruments includes instrumentation used for dental implant surgery. These instruments are reusable devices.

Intended Use

Bicon Dental Implant Surgical Instruments are used for dental implant surgery, laboratory, or restorative procedures.

Contraindications

Bicon Dental Implant Surgical Instruments should not be used for individuals with known allergic sensitivity to stainless steel or titanium aluminum vanadium alloy (Ti6Al4V) as it may cause hypersensitivity.

Warnings and Precautions

- The device is to be used on the instruction of, or by a dental or medical practitioner.
- When handling sharp instruments use extreme caution to avoid injury.
- Always inspect the device before use for worn cutting edges as bent or damaged instruments could cause device failure and result in patient or user injury.
- Proper irrigation is required when using instruments at high speed (75+ RPM); undesirable heat may be generated which can cause patient discomfort, tooth or tissue necrosis, or patient burns.
- Unless otherwise indicated, instrument sets are NOT sterile and must be sterilized prior to use.
- Do not use any other method of cleaning or sterilization, such as automated cleaning methods, other than the manual method listed in this instruction or solutions such as hydrogen peroxide, bleach, or alcohol, as this may damage the device and cause premature failure.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established, if applicable.

Residual Risks and Undesirable Effects

Residual risks and undesirable effects include adverse reaction to materials the instruments are made from, infection from a contaminated instrument, or instrument breakage during use.

Care and Handling of Instruments

Surgical instruments and instrument cases are susceptible to damage from prolonged use, and through misuse or rough handling. Care must be taken to avoid compromising their exacting performance. With proper care and handling, latch instruments are capable of being re-sterilized multiple times and may be used at least 20 times.

To minimize damage, the following should be done:

- Inspect the instrument case and instruments for damage when purchased and after each use and cleaning. Incompletely cleaned instruments should be re-cleaned, and those that need repair set aside for repair service or returned to the manufacturer or distributor.
- Only use an instrument for its intended purpose.
- When handling sharp instruments use extreme caution to avoid injury; consult with an infection control practitioner to develop and verify safety procedures appropriate for all levels of direct instrument contact.

General Instructions

Clean and sterilize non-sterile instruments in accordance with the instructions prior to first use and prior to each reuse or disposal. When inserting latch instruments into a handpiece, check to make sure that the instrument is fully seated prior to operation; damage to the instrument or to the handpiece may result from instruments not engaged fully. Usage is dependent on the surgical necessities; consult the surgical manual for proper instrument speeds if applicable. Typically, Bicon instruments, including reamers, are run at low speeds (50 RPM).

Cleaning and Decontamination

1. Disassembly – The majority of surgical instruments are simply constructed and will not require disassembly. However, some of the more complex instruments are made of several components and these should be disassembled into their individual parts prior to cleaning or decontamination.

2. Removal of Gross Contamination – The effectiveness of subsequent decontamination processes depends on prior removal of gross soil as it may be impaired by dried or coagulated protein. Clean instruments as soon as possible after use. Pre-soak instruments in an enzymatic detergent immediately after the surgical procedure for a minimum of 20 minutes. Reamers and other small parts should soak separately from the surgical kit. Ensure contact of the detergent with all instrument surfaces, taking special care with lumens, blind holes, or cannulas. Surgical instruments and kit should then be cleaned manually with a soft brush to remove all debris. To thoroughly clean instruments with multiple parts, separate into individual components. Lumens, blind holes, and cannulas should be cleaned using a snug-fitting soft brush with a twisting motion while pushing in and out multiple times. Rinse in tap water for a minimum of one minute. Thoroughly flush lumens, blind holes, and cannulas and other difficult areas.


3. Clean Ultrasonically – Clean instruments and kit in an ultrasonic machine filled with an anti-corrosive enzymatic detergent for a minimum of 10 minutes. Make sure all surfaces of the instrument are contacting the detergent. Compatible detergents and rinse aids may be used as recommended by the manufacturer of the ultrasonic unit, but these should be of neutral or near neutral pH. Excessively acidic or alkaline solutions may corrode instruments or instrument cases. Remove the instruments from the ultrasonic machine and rinse with purified water for a minimum of one minute or until there is no sign of residue detergent or debris. Make sure all areas of the instrument are rinsed.

4. Dry and Reassemble – Dry instruments with a clean, absorbent non-shedding lint-free cloth. Clean, filtered compressed air may be used to dry difficult to reach areas such as lumens, blind holes, or cannulas. Visually inspect the instruments for damage or wear that would prevent proper operation. This includes if the instrument is broken, has worn edges, shows signs of corrosion, or reference markings are illegible. If the testing is acceptable, place all instruments and small parts back into the surgical kit, if applicable. Place an autoclave indicator strip in the kit. Close the kit.

5. Sterilization – Double wrap instrument kit with autoclave paper wrap and seal with autoclave indicator tape. Individual instruments may be packaged in an approved medical grade sterilization pouch (Class 4) or wrap. Sterilize in an autoclave for a minimum of 8 minutes at 273°F (134°C) using moist heat/steam. Run the autoclave dry cycle for a minimum of 30 minutes.

Storage

Instruments cases that have been processed and wrapped to maintain sterility should be stored in a manner to avoid extremes in temperature and moisture. Care must be exercised in handling of wrapped cases to prevent damage to the sterile barrier. The health care facility should establish a shelf life for wrapped instrument cases, based upon the type of sterile wrap used and the recommendations of the sterile wrap manufacturer. The user must be aware that maintenance of sterility is event-related and that the probability of occurrence of a contaminating event increases over time and with handling, whether woven or non-woven materials, pouches, or container systems are used as the packaging method.

Key	
	Manufacturer
	Date of manufacture
	Authorized representative in the European Community
	Consult instructions for use
	By prescription only

	Batch code
	Catalogue number
	Use-by date
	Quantity
	Caution

	Do not use if package is damaged
	Do not re-use
	Non-sterile
	Do not re-sterilize
	Medical device

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