**MADE IN THE USA**

**Bicon, LLC**  
501 Arborway  
Boston, MA 02130  
USA  
tel: 617.524.4443  
[support@bicon.com](mailto:support@bicon.com)

**Bicon Europe, Ltd.**  
Unit 4 Ballycummin Village  
Ballycummin, Limerick  
Ireland  
tel: +353 6130 3494  
[ireland@bicon.com](mailto:ireland@bicon.com)

**Bicon Marketing Limited**  
1 Abacus House Newlands Road  
Corsham, Wiltshire SN13 0BH  
United Kingdom  
[uk@bicon.com](mailto:uk@bicon.com)

### **Subsidiaries\***

**Bicon Europe, Ltd.**  
Limerick, Ireland  
[ireland@bicon.com](mailto:ireland@bicon.com)

**Bicon Japan Co. Ltd.**  
Tokyo, Japan  
[japan@bicon.com](mailto:japan@bicon.com)

**Bicon Italia**  
Cisterna di Latina, Italia  
[italy@bicon.com](mailto:italy@bicon.com)

### **Australian Sponsor**

**Emergo Australia**  
201 Sussex Street • Darling Park, Tower II, Level 20  
Sydney, NSW 2000 Australia

\*For complete address information of Bicon subsidiaries and other regional distributors, please visit [www.bicon.com](http://www.bicon.com)

## General Information - Implants and Abutments

### Description

The Bicon dental implant or abutment is made from surgical grade titanium alloy (Ti6Al4V). This material has been shown to be biocompatible by many researchers. The surface configuration has been carefully designed and should not be modified. The implants are supplied grit blasted (Integra-Ti™) and hydroxylapatite coated (Integra-CP™).

Bicon implants can be used in the mandible or the maxilla. Various sizes are available to fit into created sockets in healed bone or into fresh extraction sites. For surgical procedure information, consult the surgical manual or view Bicon technique videos at <https://www.bicon.com>.

### Intended Purpose

Bicon implants and abutments are surgically placed by dentists in the jawbone to help restore the chewing function of patients who may be missing one or more teeth.

### Indications

The Bicon implant is designed for use in edentulous sites in the mandible or maxilla for support of a complete denture prosthesis, a final or intermediate abutment for fixed bridgework or for partial dentures, or as a single tooth replacement.

### Contraindications

Bicon implants should not be used in patients with contraindicated diseases such as blood dyscrasias, uncontrolled diabetes, hyperthyroidism, bruxism, oral infections, or malignancies. Nor should Bicon implants be used in patients with contraindicated conditions such as myocardial infarction within the past year or insufficient surrounding bone to permit the use of an appropriately sized implant. The 3.0 x 6.0mm, 3.0 x 8.0mm, and 3.5 x 8.0mm implants are not indicated for use as a single tooth replacement either splinted or unsplinted in the molar region. Implants should not be placed if there is insufficient alveolar bone width and height to surround the implants.

### Warnings and Precautions

Implant surgery is a highly specialized and complex procedure. Special training is required in established techniques of oral implantology. Implant courses and seminars are strongly recommended before surgery is attempted. Improper technique can result in implant failure and substantial loss of surrounding bone. Bicon implants should not be used in sites or situations other than those specifically indicated. To do so may result in failure of the implant with the concomitant destruction of supportive bone. Bicon implants are intended to be used only with the specially designed bone drills supplied by Bicon. The 3.0mm diameter implant is designed to be used with straight abutments only. When using Bicon implants, clinicians should closely monitor patients for any of the following conditions: peri-implant bone loss, changes to the implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If these implants show mobility or show greater than 50% bone loss, the implant should be evaluated for possible removal. Placing a permanent abutment on an implant that is not fully healed or is mobile may cause premature failure of the implant.

Panoramic radiographs should be used to determine if adequate bone is present at the proposed implant site as well as to locate established critical anatomical features such as the mandibular canals, mental foramina, maxillary sinuses, and adjacent teeth. Palpation and direct visual inspection of the prospective implant site are also necessary to determine the anatomy of available bone. A thorough clinical evaluation is imperative. Proper patient motivation is essential, if the procedure is to be successful. This may include having the patient stop smoking or engage in proper oral hygiene.

### MRI Safety Information

Concerning MR safety, due to the large variety of MRI scanners available on the market, Bicon cannot make any predictions regarding the safety or behavior of our implants and components in any specific MRI system. The risk assessment concluded that the devices in the Bicon Dental Implant and Abutment System can generally be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)
- Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of continuous scanning (i.e., per pulse sequence) in the Normal Operating Mode

### Residual Risks and Undesirable Effects

One hundred percent implant success cannot be guaranteed. Failure to observe the limitations of use may result in failure. Reuse of single-use devices increases risk of rejection, infection, and disease transmission. Factors such as infection, disease, inadequate bone quality and/or quantity, and patient behavior such as tobacco use or poor oral hygiene can result in osseointegration failures following surgery or initial osseointegration.

Possible complications following the insertion of dental implants include pain, swelling, bleeding, dehiscence, delayed healing, paresthesia, edema, hematoma, infection, inflammation, and generalized allergic reaction. More persistent symptoms include chronic pain in connection with implants, permanent paresthesia, nerve damage, loss of bone, infection, and fracture of the implant or prosthesis. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established, if applicable. A summary of safety and clinical performance may be found on Eudamed through the European website (<https://ec.europa.eu/tools/eudamed>) or at <https://www.bicon.com/safety>.

### Healing Plug Inserters

The Bicon healing plug inserters are made from high density polyethylene (HDPE).

### Abutments

Permanent abutments should be placed after the clinician determines whether there is sufficient bone density to support the restoration. Abutment post sizes should match the implant well size for compatibility. For example, a 2.0mm post abutment will be compatible with a 2.0mm well implant, a 2.5mm post abutment with a 2.5mm well implant, and a 3.0mm post abutment with a 3.0mm well implant. Temporary abutments, made from surgical grade titanium alloy (Ti6Al4V) or PEEK, are prosthetic components connected to the dental implant intended for use as an aid in prosthetic rehabilitation and are not intended to handle significant loads. The temporary abutment may be modified to develop a gingival sulcus or a desired shape. PEEK Temporary Abutments are intended for temporary use up to 180 days during the healing period of dental implant placement and will be replaced with a permanent abutment, such as the Universal, Stealth, Non-Shouldered, Fixed-Detachable, or Brevis abutment.

### How Supplied

Bicon dental implants and abutments are supplied in a sterile package to be stored at room temperature (15-25°C). If the package is damaged and an implant becomes contaminated, it should be discarded. The material is non-hazardous and can be disposed of in accordance with regional and local authority requirements. Bicon implants are single-use devices and should be inserted immediately after their sterile package has been opened. **Implants and abutments are provided sterile via gamma radiation and must not be reused to prevent rejection, infection, or disease transmission.**

### Sterile Handling Procedures

1. The assistant removes the sterile Tyvek blister pack from the cardboard folder and carefully opens the Tyvek blister pack and allows the inner pack to fall freely on to a sterile tray. Do not contaminate the inner pack.
2. The dentist grasps the implant component through the envelope, cuts one end of the inner pack with a pair of scissors, and removes the component from the envelope in the following manner:
  - Implant:** Grasp the healing plug inserter with gloved fingers or forceps and insert the implant into the prepared osteotomy according to Bicon's surgical instructions.
  - Abutment:** Remove the abutment from the inner pack by grasping the head portion of the abutment with a pair of forceps and insert the abutment into the fully healed well flushed implant.
3. Done properly, no one need touch sterile components directly.

### Warranty and Limitations

Bicon warranty and limitations may be found on the website at <https://www.bicon.com/warranty>.

**CAUTION:** USA Federal law restricts this device to sale by or on the order of a licensed dentist.

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