# SynthoGraft®

Pure Phase Beta-Tricalcium Phosphate

EN ENGLISH INSTRUCTIONS FOR USE

## Description

SynthoGraft is a synthetic, biocompatible and resorbable granulate ceramic made of pure-phase beta-tricalcium phosphate ( $\text{Ca}_3(\text{PO}_4)_2$ ) which is used as a matrix for bone augmentation. SynthoGraft is available in 50-500µm and 500µm-1000µm particle sizes in vials ranging from 0.25g to 2.0g. SynthoGraft is used by dentists for patients who would benefit from bone augmentation, such as edentulous or partially edentulous patients.

## Indications for Use

SynthoGraft is designed for:

- the filling and/or reconstruction of either traumatic or degenerative multi-wall bone defects.
- the augmentation of the sinus floor.
- the augmentation of alveolar ridges.
- the filling of periodontal or other alveolar bone defects and tooth sockets and osteotomies.
- the preservation of the alveolus for the preparation of an implant site.

## Product Preparation

Remove glass vial of SynthoGraft from sterile package. Dispense desired quantity of SynthoGraft particles into a sterile dish. Wet the particles with patient's blood which is most easily obtained from the ante-cubital fossa. To facilitate delivery to the surgical site, use conventional instrumentation such as a glass syringe, curette, periosteal elevator, or spatula.

## Preoperative Procedures

Dental and medical practitioners are responsible for proper preparation, surgical techniques, and post-operative care when using SynthoGraft. A proper medical history as well as a clinical and radiographic examination should be completed prior to using SynthoGraft.

## Operative Procedures

The surgical placement of SynthoGraft is equivalent to other procedures using particulate bone grafts. Conventional flap and curettage techniques are used to ensure that the site is completely debrided and that the root surfaces are thoroughly planed and smoothed. Place the previously mixed and blood-wetted SynthoGraft mixture into the surgical site using a syringe, curette, periosteal elevator, spatula or another conventional instrument. Suctioning should be limited to the periphery of the site. To reduce the possibility of particle migration, the site should not be overfilled. Placement of a resorbable collagen membrane is recommended for lateral sinus lifts only. The mucoperiosteal flaps should be sutured to achieve primary closure.

After the appropriate healing period, typically 3-6 months, a patient would be able to get a dental implant. A summary of safety and clinical performance may be found on Eudamed through the European website (<https://ec.europa.eu/>) or at <https://www.bicon.com/safety>.

## Precautions

Under no circumstances should SynthoGraft be applied in a dry state. It must be wetted with blood, since blood is necessary to ensure angiogenic closure. Similarly, SynthoGraft should not be saturated with any solution (e.g. physiological saline, NaCl, or antibiotics) other than the patient's own blood. Patients should be given an implant card allowing the identification of the device, including the device name lot number, UDI, model, and name, address, and website of the manufacturer.

## Warnings

- SynthoGraft should be used only by or under the supervision of trained personnel with experience with the surgical techniques associated with biomaterials.
- SynthoGraft is packaged and sterilized via gamma irradiation for single use only.
- Do not apply SynthoGraft dry. It must be wetted with the patient's blood.
- Do not saturate SynthoGraft with any liquid other than the patient's blood.
- Do not re-sterilize SynthoGraft.
- Discard any unused SynthoGraft granulate. Re-use of SynthoGraft may cause adverse reactions, including, but not limited to, infection, inflammation, or other injury.
- SynthoGraft is not intended for immediate loading.
- Do not overfill defect site.
- Secure SynthoGraft to prevent migration of materials.
- Do not use if package has been opened, damaged, or if the expiration date has passed.
- Do not compromise blood supply to the defect area.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established, if applicable.

## Contraindications

SynthoGraft should not be used for patients with juvenile form of chronic periodontitis, uncontrolled systemic diseases, infections, endocrinopathies, coagulation deficiencies, psychological and neurological problems, or in any other instance where the clinician believes that surgery or the use of SynthoGraft is inappropriate.

## Adverse Reactions

No adverse interactions between SynthoGraft and pharmaceuticals or other medicinal products have been reported.

## Disposal

The material is non-hazardous and can be disposed of in accordance with regional and local authority requirements.

**Rx** CAUTION: USA Federal law restricts this device to sale by or on the order of a licensed dentist.

**Bicon, LLC**  
501 Arborway  
Boston, MA 02130  
USA  
tel: 617.524.4443  
fax: 617.524.0096  
[support@bicon.com](mailto:support@bicon.com)

**Bicon Europe, Ltd.**  
Unit 4 Ballycummin Village  
Ballycummin, Limerick  
Ireland  
tel: +353 6130 3494  
fax: +353 6122-7606  
[ireland@bicon.com](mailto:ireland@bicon.com)

## Australian Sponsor

**Emergo Australia**  
201 Sussex Street • Darling Park, Tower II, Level 20  
Sydney, NSW 2000 Australia

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