# TRINIA®

The solution for **metal-free** CAD/CAM restorations

**Bicon, LLC**  
501 Arborway  
Boston, MA 02130  
USA  
tel: 617.524.4443  
[support@bicon.com](mailto:support@bicon.com)

**Bicon Europe, Ltd.**  
Unit 4 Ballycummin Village  
Ballycummin, Limerick  
Ireland  
tel: +353 6130 3494  
[ireland@bicon.com](mailto:ireland@bicon.com)

**Bicon Marketing Limited**  
1 Abacus House Newlands Road  
Corsham, Wiltshire SN13 0BH  
United Kingdom  
[uk@bicon.com](mailto:uk@bicon.com)

## Subsidiaries\*

**Bicon Europe, Ltd.**  
Limerick, Ireland  
[ireland@bicon.com](mailto:ireland@bicon.com)

**Bicon Italia, Srl.**  
Cisterna di Latina, Italia  
[italy@bicon.com](mailto:italy@bicon.com)

\*For complete address information of Bicon subsidiaries and other regional distributors, please visit [www.bicon.com](http://www.bicon.com)

## Australian Sponsor

**Emergo Australia**  
201 Sussex Street • Darling Park, Tower II, Level 20  
Sydney, NSW 2000 Australia

2797

**General**

TRINIA® CAD/CAM discs and blocks are composed of a multi-directional interlacing of fiberglass and resin in several layers. TRINIA® is milled into the desired shape of a substructure using a dental CAD/CAM system. The restorations are then bonded to tooth structure and/or to implant abutments with adhesive resin cement for use in edentulous or partially edentulous patients.

**Intended Purpose**

TRINIA® is a CAD/CAM machinable fiberglass and resin product intended to be used for the manufacture of copings, substructures, removable dentures, or frameworks.

**Indications**

TRINIA® is indicated for use in the treatment of edentulism or partial edentulism through use as copings, substructures, removable dentures, or frameworks.

**Contraindications**

TRINIA® should not be used in patients with parafunctional habits, such as bruxism.

**Warnings and Precautions**

- TRINIA® is supplied in a ready state. This material should not be fired under any circumstances.
- TRINIA® should not be milled or used at temperatures above 150°C (302°F).
- Do not contaminate TRINIA® with oils or grease.
- Do not exceed the mechanical tolerances of the device. TRINIA® has a flexural strength of >300 MPa and a flexural modulus of elasticity of <20 GPa and requires an occlusal thickness of 2.0mm for bars. Failure to observe these parameters may compromise the results achieved with TRINIA®.
- If patients are known to be allergic to any of the ingredients, TRINIA® restorations should not be used. The processing of TRINIA® discs and blocks produces dust which may irritate the skin and eyes or cause other health issues. Refer to the Material Safety Data Sheet (available at [www.bicon.com](http://www.bicon.com)).
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established, if applicable.

**Residual Risks and Undesirable Effects**

Residual risks and undesirable effects include rejection of the material or an allergic reaction to the material, or fracture of the material. A summary of safety and clinical performance may be found on Eudamed through the European website (<https://ec.europa.eu/tools/eudamed>) or at <https://www.bicon.com/safety>.

**Preparation**

Proper preparation of the tooth/teeth is essential for optimum strength, color, and retention of a milled restoration. TRINIA® is not water soluble and can be milled with air coolant by following the instructions of the manufacturer of the milling machine. Recommended bur material: Nano-diamond

**How Supplied**

TRINIA® is provided in protective packaging with a label containing its dimensions and product code. TRINIA® is provided in a ready-to-use state. Remove the product from its packaging and check for chips, cracks, or other defects.

**Storage and Use**

TRINIA® is designed to be stored in its original packaging and used at room temperature. Do not store the product in heat or in sunlight.

**Disposal**

The material is non-hazardous and can be disposed of in accordance with regional and local authority requirements.

**Rx CAUTION:** USA Federal law restricts this device to sale by or on the order of a licensed dentist or dental professional.

The solution for **metal-free** CAD/CAM restorations

**Legend**

	European Conformity		Use by date		By prescription only
	Consult instructions for use		Catalog number		Caution
	Manufacturer		Batch code		Do not use if package is damaged
	Date of manufacture		Unique device identifier		Keep away from sunlight
	Authorized representative in the European Community		Medical device		MR safe
	United Kingdom Responsible Person				

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