



DRILL STOP PILOT DRILLS

- Seven stops ranging from 5.0mm–11mm
- Self-contained in compact surgical tray
- Anodized stainless steel for increased longevity



“When the anatomy permits, Bicon’s 2.0mm diameter Drill Stop Pilot Drills allow you to confidently drill a 2.0mm pilot osteotomy confidently to a precise depth.”



2.0mm Pilot Drill Stop Kit

PART NO.	DESCRIPTION	PRICE	PART NO.	DESCRIPTION	PRICE
101-105	2.0mm Pilot Drill - 5.0mm Stop	\$80.00	101-099	2.0mm Pilot Drill Stop Kit	\$660.00
101-106	2.0mm Pilot Drill - 6.0mm Stop	\$80.00	104-999	2.0mm Pilot Drill Stop Tray	\$225.00
101-107	2.0mm Pilot Drill - 7.0mm Stop	\$80.00			
101-108	2.0mm Pilot Drill - 8.0mm Stop	\$80.00			
101-109	2.0mm Pilot Drill - 9.0mm Stop	\$80.00			
101-110	2.0mm Pilot Drill - 10mm Stop	\$80.00			
101-111	2.0mm Pilot Drill - 11mm Stop	\$80.00			